



Nursing Anne Simulator Male

Preparing future clinicians for complete patient care

Trusted by nursing educators, the Nursing Anne Simulator Male is part of the only multi-patient modular platform proven to cover every aspect of modern nursing curriculum. This high-fidelity patient simulator enables nurses to practice basic assessments, critical thinking, and advanced interventions for competence in managing gender diversity in patients and providing the highest level of care.

Available in Medium and Dark Skin Tones



SELLSHEET

www.laerdal.com

Features

Airway & Breathing

- Realistic airway with uvula, epiglottis, vocal cords and esophagus
- Various oxygen delivery methods
- Bag-mask ventilation
- Oral/nasal intubation
- Right mainstem intubation
- · Trach care and tracheal suctioning with fluid
- Spontaneous breathing synchronized with selected breath rate (0-60 bpm)
- · Compatible with ASL 5000 Lung Solution Anterior and posterior lung sound auscultation
- sites

Patient Handling

- Full articulation of arms and legs for realistic patient handling procedures
- Head can be flexed into chin to chest position and remain flexed until repositioned
- Sits unassisted and bends at waist to tripod position

Patient Assessment and Care

- Palpable anatomy for assessment and site location including clavicle, sternum, spine, anterior superior iliac spine, pubic symphysis, greater trochanter, and scapula
- · Ear canal for practice of irrigation and cleaning
- Oral care and hygiene
- Accurate anatomical landmarks for insertion of NG Tube to correct measurement
- Lavage and gavage

Medication Administration

- · Bilateral pre-ported IV arms with capability for intravenous bolus or push through intravenous infusion
- Bilateral deltoid, ventral gluteal, dorsal gluteal and thigh injection sites
- · NG tube for feeding or medication administration

Urinary & Bowel Skills/Features

- · Male genitalia with realistic anatomy
- Urinary catheterization with ability for proper fluid return
- · Enema administration and simulated rectal suppositories

Circulatory Skills/Features

- · Bilateral carotid, brachial, radial, femoral, and pedal pulses
- Pulse palpation is detected and logged
- Bilateral measurement of non-invasive blood
- pressure (auscultated or palpated) · Heart sounds synchronized with ECG

Other

- · Midclavicular central line catheter for site care and infusion procedure
- Pre-recorded vocal sounds can be played at a controllable volume
- Blinking eyes with adjustable blink rate
- Ability to open, close or partially close eyes for consciousness cue
- Interchangeable pupils (normal, dilated, constricted)
- CPR capable
- ECG monitoring capabilities when used with simulated patient monitor
- Four independently controlled auscultation areas for bowel sounds
- Complimentary clinical scenarios available on Laerdal Scenario Cloud

Patient Monitor

- · Touchscreen simulated patient monitor provides concise clinical feedback for physiological parameters
- Simulated parameters for SimPad and SimPad PLUS operation: ECG, etCO2, HR, SpO2, BP, RR, and Temperature
- LLEAP operation provides additional simulated parameters including CO2, CVP, ICP, anesthetic agent, TOF, cardiac output, and more. In addition, the LLEAP monitoring options can display x-ray, custom images, and custom video

Ordering	Information
321-05050 321-04050	Nursing Anne Simulator Male, Medium Nursing Anne Simulator Male, Dark
Operation and	Patient Monitoring Options
204-30101	SimPad PLUS (handheld remote)
400-10201	Laptop Instructor PC-Patient Monitor
400-09201	Tablet Instructor PC-Patient Monitor
400-09501	Rugged Instructor PC-Patient Monitor
400-29301	All in One Panel PC Instructor PC-Patient Monitor
Licenses	
204-50150	LLEAP for SimPad PLUS Software License
400-01050	LLEAP PC Software License
Accessories	
320-10550-M	Nursing Anne Simulator Male Kit, Medium
320-10550-D	Nursing Anne Simulator Male Kit, Dark
320-28050-M	Male Genitalia Kit; Medium
320-28050-D	Male Genitalia Kit, Dark
320-34050-M	Female Genitalia, Medium
320-34050-D	Female Genitalia, Dark
212-29650	Headset and Mic with USB (LLEAP Operation)
200-30950	SimPad Headset
185-10050	ShockLink System (adapter cables may be required)
375-51001-M	Multi-Venous Arm, Adult Female, Left
213-20050	Blue Pupil Inserts
320-27050-M	Mastectomy Module Medium
320-22550	Subcutaneous Injection Pads, 50
320-23150	Silicone Dressing, IV - 50 Pack
320-23250	Silicone Dressing, 4x4 - 50 Pack
320-23350	Silicone Dressing, 8x10 - 10 Pack

- 320-24050-M Wound Assessment and Care Kit Medium
- 320-24050-D Wound Assessment and Care Kit; Dark
- 320-25550-M Staged Wound Feet Set; Medium

320-25550-D Staged Wound Feet Dark

Take Your Simulation Training Further

Simulation based solutions and services to take competency in caring for patients further inside and outside the classroom

NextGen vSim® for Nursing

Adaptive, interactive virtual simulation learning experience



ASL 5000[™] Lung Solution Highly realistic ventilation management training



Laerdal Scenario Cloud

Digital library with a variety of expert-validated clinical scenarios



SimCapture Simulation and Debriefing Management Solution



Operation

Nursing Anne Simulator Male can be operated by a SimPad®/SimPad® PLUS remote or a LLEAP instructor computer, can be used with a simulated patient monitor, and is compatible with an ASL 5000 Lung Solution.





For more information, visit www.laerdal.com